



iPrep[™] PureLink[™] Virus Kit

**For purification of viral nucleic acids from
cell-free samples using the iPrep[™]
Purification Instrument**

Catalog no. IS10008

Rev. date: 15 August 2010
Manual part no. 100004281

MAN0000409

User Manual

Contents

Experienced Users Procedure	iv
Kit Contents and Storage	vi
Introduction	1
Product Overview	1
iPrep™ Purification Instrument	3
Methods.....	6
General Information	6
Isolating RNA/DNA from Cell-free Samples	8
Expected Results.....	14
Troubleshooting	17
Appendix	19
Accessory Products	19
Technical Support	20
Purchaser Notification.....	22

Experienced Users Procedure

Introduction This quick reference sheet is included for experienced users of the iPrep™ PureLink™ Virus Kit. For more details, refer to this manual.

Step	Procedure
Purification Protocol	<div><div>1. Mix fresh cell-free samples or thaw frozen cell-free samples or lyse bacteria using lysozyme for bacterial gDNA isolation (page 9).</div><div>Note: The need for lysozyme digestion step is dependent on the type of bacteria in use as certain Gram negative bacteria (e.g., <i>E. coli</i>) do not need the lysozyme digestion step. If lysozyme digestion step is not performed, use 200 µL or 400 µL bacterial samples.</div><div>2. Open the iPrep™ Card Slot and insert the iPrep™ Card: Viral DNA/RNA in the slot (arrow on the card is at the top and card label is facing your left side).</div><div>3. Turn ON the iPrep™ Instrument using the power switch on the left side of the instrument.</div><div>The digital display shows the version for the iPrep™ which changes in few seconds to display the Main menu.</div><div>4. Press Start to run a protocol.</div><div>5. Open the iPrep™ instrument door and remove iPrep™ Racks to set up the platform.</div><div>6. Remove the iPrep™ PureLink™ Virus Cartridges from the box. To collect any solution from the foil, tap the cartridge to deposit the solution at the bottom of the tube.</div><div>7. Insert one iPrep™ Sample Processing Tube (2 mL) in the heated tube position of cartridge (position 11) for each iPrep™ PureLink™ Virus Cartridge that is used.</div><div>8. Load the cartridges on the iPrep™ Cartridge Rack and insert the loaded rack on to the iPrep™ Platform.</div></div>

Continued on next page

Experienced Users Procedure, Continued

Step	Procedure
Purification Protocol, Continued	<ol style="list-style-type: none"> 9. Load the iPrep™ Tip and Tube Rack as follows: <ul style="list-style-type: none"> • Load the first row (labeled as E) with 1–13 elution tubes without caps • Load the second row (labeled as T1) with iPrep™ Small Tips (blue tips, page 19) in iPrep™ Tip Holders for eluting samples in 20 µL elution volume • Load the third row (labeled as T2) with iPrep™ Tips in iPrep™ Tip Holders • Load the fourth row (labeled as S) with iPrep™ Sample and Elution Tubes containing samples 10. Read the sample and elution tube barcodes, if needed. 11. Insert the iPrep™ Tip and Tube Rack on to the iPrep™ Platform. 12. Close the iPrep™ instrument door. Press Enter (↵) to continue. 13. When prompted, select the appropriate lysis mode, sample volume, and elution volume. 14. Press Start. The automated purification protocol begins and various steps of the protocol including the approximate time remaining are displayed on the digital display. 15. At the end of the run, the instrument beeps briefly and the digital display shows Protocol Finished for 10 seconds. The Main menu appears after 10 seconds. 16. Open the instrument door. 17. Remove and cap the elution tubes containing the purified nucleic acid. Store the purified RNA/DNA at –80°C. 18. Discard the used cartridges, tips, and tubes into biohazard waste. Do not reuse the cartridges. 19. To purify more samples using the same iPrep™ Card, load the racks with new cartridges, tips, tubes, and samples, and start the protocol as described. 20. If you are not using the instrument, close the instrument door and turn the power switch to OFF. 21. Remove the iPrep™ Card and store card in the box.

Kit Contents and Storage

Shipping and Storage

The iPrep™ PureLink™ Virus Kit is shipped at room temperature.

Upon receipt, store the iPrep™ PureLink™ Virus Kit at room temperature. See below for kit contents.

All components are guaranteed stable for 6 months when stored properly.

Kit Contents

The components supplied in the iPrep™ PureLink™ Virus Kit are listed below.

Sufficient reagents are supplied to perform 52 purifications.

Reagents	Amount
iPrep™ PureLink™ Virus Cartridge Kit	4 cartridge racks
iPrep™ Sample and Elution Tubes	2 bags with 52 tubes
iPrep™ Sample Processing Tubes	1 bag with 52 tubes
iPrep™ Tips and Tip Holders	1 bag with 52 tips and holders

Continued on next page

Kit Contents and Storage, Continued

**iPrep™
PureLink™
Virus
Cartridge Kit
Contents**

Each iPrep™ PureLink™ Virus Cartridge has 12 positions with 10 sealed wells and two heating positions (position 12 with an empty well and position 11 to add an empty or reagent filled tube).

Positions 1–10 contain wells filled with reagents for this protocol.

The components supplied in each well of the iPrep™ PureLink™ Virus Cartridge Kit are listed below.

Store the iPrep™ PureLink™ Virus Cartridge Kit at room temperature and do not freeze the cartridge kit.

Reagent	Well no.
Viral Lysis Buffer	1
Proteinase K (20 mg/mL) in storage buffer (proprietary)	2
Elution Buffer (DEPC Water)	3
Empty	4
Wash Buffer 2	5
Wash Buffer 2	6
Wash Buffer 1	7
Dynabeads® MyOne™ SILANE (4.8 mg/mL in Bead Storage Buffer)	8
Rinse Buffer	9
Isopropanol	10

Intended Use

For research use only. Not intended for any animal or human therapeutic or diagnostic use.

Introduction

Product Overview

Introduction

The iPrep™ PureLink™ Virus Kit allows rapid and automated extraction of viral RNA/DNA as well as bacterial genomic DNA (gDNA) from fresh or frozen cell-free biological fluids (e.g., plasma, serum, urine) and cell culture supernatants.

Nucleic acid is purified from viral or bacterial cell-free samples using the Dynabeads® MyOne™ SILANE and iPrep™ Purification Instrument within 40 minutes without the use of centrifugation.

The purified RNA/DNA is devoid of proteins and nucleases, and is suitable for use in downstream applications that allow viral detection and genotyping.

iPrep™ Purification Instrument

The iPrep™ PureLink™ Virus Kit is designed for use with the iPrep™ Purification Instrument.

The iPrep™ Purification Instrument is a benchtop, automated nucleic acid purification instrument with integrated Magnetic and Syringe Unit capable of purifying nucleic acids from up to 13 samples (12 samples + 1 positive control) using magnetic bead-based technology. See page 3 for details on the iPrep™ Purification Instrument.

System Overview

The iPrep™ PureLink™ Virus Kit combines the sensitivity and capacity of Dynabeads® MyOne™ SILANE with the speed and convenience of the iPrep™ Instrument to allow automated purification of high-quality DNA and RNA from up to 13 samples (12 samples + 1 positive control) within 40 minutes. The Dynabeads® MyOne™ SILANE are monodisperse magnetic beads (1 µm) with an optimized silica-like surface chemistry and a high specific surface area. Purification is achieved using magnetic bead-based procedure, and avoids the use centrifuges or vacuum manifolds.

The viral particles are lysed using Viral Lysis Buffer and proteins are digested with Proteinase K. The lysate is mixed with Dynabeads® MyOne™ SILANE for subsequent nucleic acid binding to the beads. The nucleic acid-bound magnetic beads are separated from the lysate using magnetic separation. The beads are thoroughly washed with Wash Buffers to remove contaminants. The RNA/DNA is eluted in Elution Buffer.

Continued on next page

Product Product Overview, Continued

Advantages

The iPrep™ PureLink™ Virus Kit provides the following advantages:

- Uses a magnetic bead-based technology to isolate viral RNA/DNA and bacterial gDNA without the need for centrifugation or vacuum manifolds
- Rapid and automated purification of nucleic acids within 40 minutes from a wide range of viral samples as well as Gram positive and Gram negative bacteria using the iPrep™ Instrument
- Pre-filled reagent cartridges provide easy set up and consistent results
- Minimal sample cross-contamination
- Purified nucleic acids demonstrate improved downstream performance in various applications

System Specifications

Starting Material:	200 µL or 400 µL cell-free sample (see page 6)
Bead Size:	~1 µm
Bead Amount per Reaction:	2.4 mg
Number of Samples:	Up to 13
Elution Volume:	20 µL*, 50 µL, or 100 µL

*To elute samples using 20 µL elution volume, you need to order iPrep™ Small Tips (page 19).

iPrep™ Purification Instrument

Introduction

The iPrep™ Purification Instrument is a benchtop, automated nucleic acid purification instrument with integrated Magnetic and Syringe Unit capable of purifying nucleic acids from up to 12 samples and one positive control. Each iPrep™ Instrument consists of the Magnetic and Syringe Unit, and a platform. A pre-programmed iPrep™ Protocol Card controls the purification parameters such as buffer volumes, mixing steps, and incubation time. For more details on the iPrep™ Purification Instrument, see the manual supplied with the instrument.

iPrep™ Reaction Cartridge


The iPrep™ Reaction Cartridges are supplied with iPrep™ Kits and are designed to fit onto the iPrep™ Cartridge Rack in only one orientation. Each cartridge is pre-filled with reagents required for the iPrep™ PureLink™ Virus protocol.

Each cartridge has 12 positions with 10 sealed wells and two heating positions (position 12 with an empty well and position 11 to add an empty or reagent filled tube). For the iPrep™ PureLink™ Virus Kit, positions 1–10 contain wells filled with reagents.

Cartridge Specifications:

Material:	Polypropylene cartridge sealed with laminated aluminum foil
Max Volume:	1,000 µL/well
Dimension:	5.9 inches (l) × 1.2 inches (w) × 0.7 inches (d)

Note: The image below shows an example of an iPrep™ Reagent Cartridge and is not an image of an iPrep™ Viral Cartridge.



Continued on next page

iPrep™ Purification Instrument, Continued

iPrep™ Tips and Tip Holders

The iPrep™ Tips and Tip Holders are included with iPrep™ Kits and are placed on the iPrep™ Tip and Tube Rack as described on page 11. While assembling tips on the rack, insert the iPrep™ Tips into the iPrep™ Tip Holders using gloved hands. Always use tips with the holders to prevent any contamination.

iPrep™ Small Tips are available separately from Invitrogen (page 19) and are used for lower elution volumes.


Tip Specifications:

Tip Material:	Polypropylene with filter barriers
Tip Holder Material:	Polypropylene
Volume:	5–1,000 µL (standard tips) 2–200 µL (small tips)
Standard Tip Dimensions:	3.9 inches (l) × 0.43 inches (d)
Small Tip Dimensions:	2.9 inches (l) × 0.43 inches (d)

iPrep™ Tip Holder



iPrep™ Tip



Continued on next page

iPrep™ Purification Instrument, Continued

iPrep™ Tubes

Two sets of iPrep™ Tubes are required for the purification protocol. The iPrep™ Sample and Elution Tubes and iPrep™ Sample Processing tubes are included with each iPrep™ Viral Kit and placed on the iPrep™ Tip and Tube Rack as described on page 11.

Tube Specifications:

Material:	Polypropylene
Capacity:	1.5 mL (iPrep™ Sample and Elution Tubes) 2.0 mL (iPrep™ Sample Processing Tubes)
Style:	Tubes with caps
Dimensions:	1.7 inches (l) × 0.4 inches (d)



iPrep™ Card: Virus

To isolate viral RNA/DNA and bacterial gDNA using the iPrep™ PureLink™ Virus Kit with the iPrep™ Purification Instrument, you need to purchase the iPrep™ Card: Viral DNA/RNA (page 19).

The iPrep™ Card: Viral DNA/RNA is pre-programmed with the purification protocol for cell-free samples that directs the volume of reagents used and incubation time.

Always store the card in the box, protected from light.

To avoid damaging the card:

- Do not drop or bend the card
- Do not wipe or clean the card using volatile chemicals such as alcohol or equivalent
- Do not expose the card to water

iPrep™ Platform

The platform on the iPrep™ Instrument allows the placement of iPrep™ Tip and Tube Rack, and iPrep™ Cartridge Rack that are filled with plastic disposables and reagent cartridges required for the purification protocol.

Set up the platform as shown in the figure on page 11 for the iPrep™ PureLink™ Virus Kit.

Methods

General Information

User Supplied Materials

In addition to the reagents supplied with the kit, you also need the following materials and instrumentation:

- iPrep™ Purification Instrument (page 19)
 - iPrep™ Card: Viral DNA/RNA (page 19)
 - Viral or bacterial samples (see below)
-

Cell-free Samples

The iPrep™ PureLink™ Viral Kit is designed to isolate viral and bacterial nucleic acids from cell-free biological fluids such as plasma, serum, and urine as well as cell culture supernatants using fresh or frozen samples.

To obtain high yield of nucleic acids and minimize any degradation, follow these guidelines:

- Collect the sample (such as plasma or serum) and proceed immediately to the purification protocol (page 10). If desired, you can store the sample at 4°C for short-term storage (up to 4 hours) or freeze the sample at -20°C or -80°C for long-term storage.
 - Do not freeze-thaw the plasma or serum sample more than once.
 - Remove any visible cryoprecipitates from samples by centrifugation at $\sim 7,000\times g$ for 2-3 minutes. Use the clear supernatant immediately for purification.
 - If you need to concentrate the cell culture supernatant use appropriate centrifugal concentrators.
-

Sample Volume

The iPrep™ PureLink™ Viral Protocol is designed to purify nucleic acids from **200 μL or 400 μL** cell-free samples.

Do not use less than 200 μL sample volume as using less sample volume results in excessive bubble formation during the purification protocol thereby lowering the yield.

Continued on next page

General Information, Continued

Safety Information

Follow the safety guidelines below when using the iPrep™ PureLink™ Virus Kit.

- Treat all reagents supplied in the kit as potential irritants.
 - Always wear a suitable lab coat, disposable gloves, and protective goggles when handling samples.
 - Dispose of viral and bacterial cell-free samples as biohazardous waste.
-



Follow the recommendations below to obtain the best results:

- Use disposable, individually wrapped, sterile plastic ware
 - Use only sterile, new pipette tips (aerosol-barrier pipette tips recommended) and microcentrifuge tubes
 - Wear latex gloves while handling reagents and RNA samples to prevent RNase contamination from the surface of the skin
 - Always use proper microbiological aseptic techniques when working with RNA
 - Use RNase AWAY® Reagent (page 19) to remove RNase contamination from surfaces
 - **Do not** freeze the beads as this irreparably damages them. Store the beads at room temperature.
 - When using beads from the Reaction Cartridges, collect any solution from the foil by tapping the cartridge to deposit the solution at the bottom of the tube. **Do not** allow the beads to dry out as this renders them non-functional.
 - Discard Reaction Cartridges, iPrep™ Tips, and iPrep™ Tip Holders after use. **Do not reuse.**
-

Isolating RNA/DNA from Cell-free Samples

Introduction

Instructions to isolate viral RNA/DNA and bacterial genomic DNA from cell-free samples using the iPrep™ PureLink™ Virus Kit with the iPrep™ Purification Instrument are described below.

Starting Material

Use this procedure to isolate viral RNA/DNA and bacterial genomic DNA from **200 µL or 400 µL** cell-free samples. See page 6 for sample volume.

Materials Needed

- Cell-free samples (such as plasma or serum samples, page 6)
- iPrep™ Purification Instrument (page 19)
- iPrep™ Card: Viral DNA/RNA (page 19)
- *Optional:* iPrep™ Small Tips (page 19)
- *Optional:* Carrier RNA (page 9)
- *Optional:* Lysozyme Digestion Buffer (25 mM Tris-HCl, pH 8.0, 2.5 mM EDTA, 1% Triton® X-100) and lysozyme for bacterial gDNA isolation

Components Supplied with the Kit

- iPrep™ PureLink™ Virus Cartridge Kit
 - iPrep™ Sample and Elution Tubes
 - iPrep™ Sample Processing Tubes
 - iPrep™ Tips and iPrep™ Tip Holders
-

Before Starting

Perform the following before starting:

- Thaw frozen cell-free samples or mix the fresh cell-free samples or lyse bacteria using lysozyme for bacterial gDNA isolation (page 9), and store on ice until use
 - Ensure that you have the iPrep™ Card: Viral DNA/RNA (page 19) to run the protocol
 - Make sure the iPrep™ Purification Instrument is unpacked and installed
-

Continued on next page

Isolating RNA/DNA from Cell-free Samples, Continued

Elution Volume

The iPrep™ PureLink™ Viral Kit utilizes low, recommended elution volume of 50 µL and 100 µL to elute highly concentrated nucleic acids required for sensitive downstream applications.

You may elute the nucleic acids in an elution volume of 20 µL depending on your downstream application.

Note: You need to purchase iPrep™ Small Tips available separately from Invitrogen (page 19) to use 20 µL elution volume.

Carrier RNA

The iPrep™ PureLink™ Viral Kit protocol does not include carrier RNA during the lysis step as we have not observed any detectable improvement in the performance with carrier RNA.

However, if you wish to include carrier RNA during lysis, use 5.6 µg of carrier RNA (see below for **Sample Lysis**).

Sample Lysis

The iPrep™ PureLink™ Viral Kit protocol includes an **inline** lysis step, i.e., the sample is placed in the tube on the rack and once the protocol is started, sample lysis occurs in the iPrep™ whereby the Viral Lysis Buffer and Proteinase K are added to the sample and lysis performed at 65°C for 10 minutes. The inline lysis step reduces sample handling and exposure to any infectious material.

However, if you wish to perform **offline** lysis, i.e., lyse the sample prior to placing samples in the iPrep™, perform lysis using the PureLink™ Viral Lysis Buffer and Proteinase K available separately from Invitrogen (page 19). For best results, use only PureLink™ Viral Lysis Buffer and Proteinase K and not any other lysis buffer for sample lysis.

If you wish to use carrier RNA, use 5.6 µg carrier RNA (Yeast tRNA, available from Invitrogen, page 19) to the samples. Perform offline lysis of samples in the presence of carrier RNA or load the samples with carrier RNA to the iPrep™ Instrument for inline lysis.

Note: If you perform offline lysis, the instrument automatically advances to the next step in the protocol and only performs the lysis incubation step.

Continued on next page

Isolating RNA/DNA from Cell-free Samples, Continued

Preparing Bacterial Lysates

To isolate bacterial gDNA from cell-free samples, prepare the lysates as follows and load the samples to the iPrep™ Instrument.

Note: The need for lysozyme digestion step is dependent on the type of bacteria in use as certain Gram negative bacteria (e.g., *E. coli*) do not need the lysozyme digestion step. If lysozyme digestion step is not performed, load 200 µL or 400 µL bacterial samples to the iPrep™ Instrument.

1. Prepare Lysozyme Digestion Buffer (25 mM Tris-HCl, pH 8.0, 2.5 mM EDTA, 1% Triton® X-100).
2. To 200 µL of Lysozyme Digestion Buffer, add fresh lysozyme to a final concentration of 20 mg/mL.
3. To 200 µL of cell-free sample, add 200 µL 20 mg/mL lysozyme in Lysozyme Digestion Buffer from Step 2. Mix thoroughly.
4. Incubate samples at 37°C for 30 minutes.
5. Mix briefly again following incubation and add 400 µL sample to the iPrep Sample and Elution Tube (page 11).

Purification Protocol

Purify viral RNA/DNA and bacterial genomic DNA from cell-free samples using the iPrep™ Purification Instrument as described below.

For details on using the iPrep™ Purification Instrument, refer to the manual supplied with the instrument.

Insert the iPrep™ Card: Viral DNA/RNA (available separately from Invitrogen, page 19) prior to turning on the instrument.

1. Ensure the power switch on the iPrep™ Instrument is on the OFF position.
2. Open the iPrep™ Card Slot and insert the iPrep™ Card into the slot in the correct orientation (arrow on the card is at the top and card label is facing your left side).
3. Using the power switch located on the left side of the instrument, turn ON the instrument.

If the card is fully inserted in the correct orientation, all axes return to their original positions automatically. The digital display shows the version for the iPrep™ which changes in a few seconds to display the Main menu.

Continued on next page

Isolating RNA/DNA from Cell-free Samples, Continued


Purification Protocol, Continued

Procedure continued from previous page

4. Press **Start** to run a protocol.
5. Open the iPrep™ instrument door. Remove the iPrep™ Cartridge Rack, and iPrep™ Tip and Tube Rack to set up the platform.
6. Remove the desired number of iPrep™ PureLink™ Virus Cartridges from the box. To collect any solution from the foil, tap the cartridge to deposit the solution at the bottom of the tube.

Note: You can load 1–13 cartridges on the rack depending on the number of samples that you wish to process. If you are loading less than 13 cartridges, ensure that the remaining plastic ware (tips and tubes) are also loaded in the same order as the cartridges.

7. Insert the iPrep™ Sample Processing Tube **in the heated tube position of the cartridge (position 11)** for each of the iPrep™ Virus Cartridge that is used.
8. Load the cartridges on the iPrep™ Cartridge Rack and insert the loaded rack on the iPrep™ platform.
9. Load the iPrep™ Tip and Tube Rack as follows:
 - Load the first row (labeled as **E**) with 1–13 elution tubes **without caps** (you may place the caps on the rack as shown in the figure below)
 - Load the second row (labeled as **T1**) with iPrep™ Small Tips (page 19) in iPrep™ Tip Holders, if you wish to use 20 µL elution volume
 - Load the third row (labeled as **T2**) with iPrep™ Tips in iPrep™ Tip Holders
 - Load the fourth row (labeled as **S**) with iPrep™ Sample and Elution Tubes **containing samples**.



Continued on next page

Isolating RNA/DNA from Cell-free Samples, Continued

Purification Protocol, Continued

Procedure continued from previous page

10. Read the sample and elution tube barcode, if needed.
11. Insert the iPrep Tip and Tube rack on the iPrep™ platform as shown on the next page.
12. Close the iPrep™ instrument door.
13. Press **Enter** (↵) to continue.
14. When prompted, select the appropriate lysis mode (inline or offline), sample volume (200 µL or 400 µL), and elution volume (20 µL, 50 µL, or 100 µL).
15. Ensure that you have loaded the cartridges, tubes, and tips in the appropriate positions, and elution tubes **do not have any caps. Make sure you have loaded a 2 mL tube in the heated tube position of the cartridge (position 11).**
16. Press **Start**. The automated purification protocol begins and various steps of the protocol including the approximate time remaining are displayed on the digital display.

Important: Do not open the door once the protocol has begun.

To pause the protocol, press the **Stop** key. To resume the protocol after a pause, press the **Start** key. To cancel/stop the protocol, press the **Stop** key twice. For details, see the iPrep™ Instrument manual.

The run times are approximate and depend on the parameters chosen such as lysis mode, sample volume, and elution volume as well as the time required to preheat the heat block.
17. At the end of the run, the instrument beeps briefly and digital display shows **Protocol Finished** for 10 seconds. The Main menu appears after 10 seconds.
18. Open the instrument door. Remove and cap the elution tubes containing the purified nucleic acid. Use the RNA/DNA for the desired downstream application or store the purified RNA/DNA at –80°C.
19. Discard the used cartridges, tips, and sample tubes into biohazard waste. **Do not reuse the cartridges.**

Continued on next page

Isolating RNA/DNA from Cell-free Samples, Continued

Purification Protocol, Continued

20. To purify more samples using the same iPrep™ Card, load the racks with new cartridges, tips, and samples, and start the protocol as described above.
 21. If you are not using the instrument, close the instrument door and turn the power switch to **OFF**.
 22. Remove the iPrep™ Card and store the card in the box, protected from light.
-

Analyzing RNA/DNA

Since the amount of viral RNA/DNA and bacterial gDNA present in cell-free samples is low we recommend that you do not use UV absorbance to determine yield. Use qRT-PCR or RT-PCR for RNA virus, and qPCR and PCR for DNA virus and bacterial gDNA using appropriate probes to determine yield or presence of nucleic acids.

To analyze viral nucleic acid size, use agarose gel electrophoresis followed by hybridization using specific labeled probes and autoradiography.

Expected Results


Adenovirus DNA

Duplicate human plasma samples (400 μ L) were spiked with 10^1 , 10^2 , 10^3 , 10^4 , and 10^5 adenoviral transfection units (TU). Viral DNA was purified using the iPrep™ PureLink™ Viral Kit as described in this manual.


qPCR was performed on the purified DNA (5 μ L) using the Platinum® Quantitative PCR SuperMix-UDG w/ROX™ dye (page 19) with TaqMan® primers in an ABI 7900 HT instrument.

Results: Detection was linear over a 5-fold range, with sensitivity down to 10 TU per sample.

Amplification Plot



Linearity Plot



Continued on next page

Expected Results, Continued


Lentivirus RNA

Duplicate human serum samples (400 μ L) were spiked with 10^1 , 10^2 , 10^3 , 10^4 , and 10^5 lentivirus transfection units (TU). Viral RNA was purified using the iPrep™ PureLink™ Viral Kit as described in this manual.


qRT-PCR was performed on the purified RNA(5 μ L) using the SuperScript™ III Platinum® One-Step qRT-PCR Kit w/ROX™ dye (page 19) with TaqMan® primers in an ABI 7900 HT instrument.

Results: Detection was linear over a 4-fold range, with sensitivity down to 10 TU per sample.

Amplification Plot



Linearity Plot



Continued on next page

Expected Results, Continued


Bacterial gDNA Results

Duplicate human plasma samples (200 μ L) were spiked with 10^2 , 10^3 , 10^4 , and 10^5 *B. subtilis* colony forming units (CFU). Bacterial gDNA was purified using the iPrep™ PureLink™ Viral Kit as described in this manual using the lysozyme digestion step.


qPCR was performed on the purified DNA (5 μ L) using the Platinum® Quantitative PCR SuperMix-UDG w/ROX™ dye (page 19) in an ABI 7900 HT instrument.

Results: Detection was linear over a 4-fold range, with sensitivity down to 100 CFU per sample.

Amplification Plot



Linearity Plot



Troubleshooting

Introduction

Refer to the table below to troubleshoot problems with the kit. To troubleshoot problems with the iPrep™ Purification Instrument, refer to the instrument manual.

Observation	Cause	Solution
Low nucleic acid yield	Too much starting material used	The purification protocol is designed for use with 200 µL and 400 µL cell-free sample volume. Using greater than the recommended amount of starting material may overload the system and cause clumping which reduces the yield.
	Incomplete lysis	Be sure to perform lysis of bacterial samples as described on page 10 prior to loading samples on the iPrep™ Instrument. If you are performing offline lysis, be sure to use PureLink™ Viral Lysis Buffer and Proteinase K (page 19) for best results.
	Insufficient amount of Dynabeads® MyOne™ SILANE added	During shipping, some Dynabeads® MyOne™ SILANE solution may adhere to the sealing foil of the cartridge. To collect any bead solution from the foil, tap the cartridge to deposit the bead solution at the bottom of the tube.
	Poor quality of sample material	Avoid repeated freezing and thawing of samples. Use fresh samples and process immediately after collection or use samples thawed only once for best results. Check the quality of the RNA in the original samples using qRT-PCR or RT-PCR.
	RNA/DNA quantitation performed using UV absorbance	Since viral nucleic acids are present in low amounts in cell-free samples, do not use UV absorbance for quantitation. Analyze viral nucleic acids using qRT-PCR, RT-PCR, qPCR, or PCR.
No nucleic acid recovered	Magnetic beads stored or handled improperly	<ul style="list-style-type: none">• Store cartridge containing the beads at room temperature. Do not freeze the cartridge as the beads may be irreparably damaged.• Make sure that the beads are in solution at all times and do not dry. Dried beads are non-functional.

Troubleshooting, Continued

Observation	Cause	Solution
Eluate is discolored	Magnetic beads present in the eluate	Remove any magnetic beads using a magnetic separator (MagnaRack™ separator is available from Invitrogen, see page 19) or centrifuge the sample in a microcentrifuge for 1 minute at maximum speed.
DNA or RNA is sheared or degraded	Bubbles formed during mixing steps	Make sure that the sample volume is at least 200 µL to prevent excessive bubble formation during mixing.
	Purified DNA repeatedly frozen and thawed	Aliquot purified DNA and store at 4°C (short-term) or –20°C (long-term). Avoid repeated freezing and thawing.
	DNA contaminated with DNases	Maintain a sterile environment while working (<i>i.e.</i> wear gloves and use DNase-free reagents).
	RNA contaminated with RNase	Follow the guidelines on page 7 to prevent RNase contamination.
Poor performance of nucleic acids in downstream enzymatic reactions	Reagents for enzymatic reactions inactive	Ensure that the enzymes and reagents used for performing downstream applications have not expired or inactivated. Repeat the reaction with fresh enzyme and reagents.
	Viral nucleic acid eluate too dilute	Optimize the amount of viral nucleic acid eluate required for your specific application and perform elution using the desired elution volume (20 µL, 50 µL, and 100 µL).

Appendix

Accessory Products

Additional Products

The table below lists additional products available from Invitrogen for use with the iPrep™ PureLink™ Virus Kit.

For more information, visit www.invitrogen.com or contact Technical Support (page 20).

Product	Amount	Cat. no.
iPrep™ Purification Instrument	1 unit	IS10000
iPrep™ Card: Viral DNA/RNA	1 card	IS10016
iPrep™ Card: gDNA Blood	1 card	IS10012
iPrep™ Card: Buffy Coat	1 card	IS10015
iPrep™ Card: gDNA Tissue	1 card	IS10013
iPrep™ Card: gDNA Forensic (includes buccal protocol)	1 card	IS10011
iPrep™ ChargeSwitch® Forensic Kit	1 kit (52 purifications)	IS10002
iPrep™ ChargeSwitch® Buccal Cell Kit	1 kit (52 purifications)	IS10003
iPrep™ ChargeSwitch® gDNA Tissue Kit	1 kit (52 purifications)	IS10004
iPrep™ Small Tips	1 bag of 52 tips	IS10111
iPrep™ Tip and Tube Rack	1 rack	IS10101
iPrep™ Cartridge Rack	1 rack	IS10102
PureLink™ Viral Lysis Buffer	500 mL	12282-500
Proteinase K	5 mL	25530-049
Yeast tRNA	25 mg	15401-011
Quant-iT™ PicoGreen® dsDNA Assay Kit	1 kit	P7589
Quant-iT™ DNA Assay Kit, High Sensitivity	1000 assays	Q33120
Quant-iT™ DNA Assay Kit, Broad-Range	1000 assays	Q33130
Qubit® Fluorometer	1 each	Q32857
MagnaRack™ Magnetic Separator	1 rack	CS15000
Platinum® Quantitative PCR SuperMix-UDG w/ROX	100 reactions	11743-100
SuperScript® III One-Step RT-PCR System with Platinum® Taq DNA Polymerase	100 reactions	12574-026
RNase AWAY™ Reagent	250 ml	10328-011

Technical Support

World Wide Web



Visit the Invitrogen website at www.invitrogen.com for:

- Technical resources, including manuals, vector maps and sequences, application notes, SDSs, FAQs, formulations, citations, handbooks, etc.
 - Complete technical support contact information
 - Access to the Invitrogen Online Catalog
 - Additional product information and special offers
-

Contact Us

For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our Web page (www.invitrogen.com).

Corporate Headquarters:

5791 Van Allen Way

Carlsbad, CA 92008 USA

Tel: 1 760 603 7200

Tel (Toll Free): 1 800 955 6288

Fax: 1 760 602 6500

E-mail:

tech_support@invitrogen.com

European Headquarters:

Inchinnan Business Park

3 Fountain Drive

Paisley PA4 9RF, UK

Tel: +44 (0) 141 814 6100

Tech Fax: +44 (0) 141 814 6117

E-mail:

eurotech@invitrogen.com

Continued on next page

Technical Support, Continued

SDS Requests

SDSs (Safety Data Sheets) are available on our website at www.invitrogen.com/sds.

Limited Warranty

Invitrogen is committed to providing our customers with high-quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, contact our Technical Support Representatives.

Invitrogen warrants that all of its products will perform according to specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order.

Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Support Representatives.

Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.

Purchaser Notification

Limited Use Label License No. 5: Invitrogen Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. For products that are subject to multiple limited use label licenses, the most restrictive terms apply. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned or controlled by Invitrogen Corporation which cover this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 5791 Van Allen Way, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500. Email: outlicensing@invitrogen.com

©2010 Life Technologies Corporation. All rights reserved.

The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners. TaqMan[®] is a registered trademark of Roche Molecular Systems, Inc. RNase AWAY is a trademark of Molecular Bio-Products, Inc. Triton is a trademark of Union Carbide Corporation.

Notes



Corporate Headquarters

5791 Van Allen Way
Carlsbad, CA 92008

T: 1 760 603 7200

F: 1 760 602 6500

E: tech_support@invitrogen.com

For country-specific contact information visit our web site at www.invitrogen.com