

Chemistry 3A

Introductory General Chemistry

Experiment 0a

Lab Exploration

Special Guidance Today

- Because you did not have access to your laboratory manuals for this start of lab session, I will quickly guide you to what you need to pay attention to
- You will have your manuals today for reference
- You have a data report form you fill out today
- Use ink pen. NO PENCILS. Just cross-out mistakes
- There are 3 parts to Expt 0a
- We skip Part 1 – Scavenger Hunt
- It has you looking for 20 items; we can learn this later

Equipment You Will Use



electronic balance



**50 mL
volumetric
flask**



**50 mL
beaker**


Equipment You Will Use

25 mL
buret

0.1 mL
division marks



100 mL
graduated cylinder



Equipment You Will Use



Part 2: Using A Buret

Equipment/Material Needed

Buret *handle carefully. It is an easily breakable analytical device!*

Buret clamp

50 mL beaker

50 mL graduated cylinder

Erlenmeyer flask

Consumables Needed

DI water (~60 mL)

food coloring (2 dyes)

- ❖ Page 12: mix 1 gtt food dye in ~30 mL DI H₂O in beaker; fill buret; record initial reading on buret; record vol + sketch
- ❖ p 13: fill grad cylinder ~20-30 mL DI H₂O; record vol + sketch; pour in Erlenmeyer; add different food dye color; mix
- ❖ from buret, add ~ 5 mL dyed H₂O to Erlenmeyer flask while swirling to mix (DON'T LET FLASK HIT BURET TIP)
- ❖ Buret stopcock control skill: 1 mL volume + droplet formation

Part 3: Using Pipets, Volumetric Flasks

Equipment/Material Needed

50 mL volumetric flask (2) w/ stopper	50 mL beaker
scoopula	10.00 mL volumetric pipet
balance	

Consumables Needed

4-6 g NaCl (sodium chloride)	food coloring (2 dyes)
DI water	

- ❖ Page 14: tare balance; weigh empty 50 mL stoppered vol. flask; record to data report
- ❖ Tare 50 mL beaker; add NaCl; record mass
- ❖ Add ~10 mL DI water; swirl to dissolve as much as possible; transfer to 50 mL vol flask; add gtt dye; add more H₂O, follow transfer method; fill to flask mark w/ transfer pipet
- ❖ Take 10.00 mL from flask to new 50 ml flask (use bulb carefully). Fill to flask mark as before using transfer pipet
- ❖ Compare color of solutions, fill out data report

Clean Up

- Water and salt solutions can safely be poured down the drain. (Like any household solutions)
This will not be the case for many solutions which may be strong acids & bases or certain inorganic compounds/salts
- Glassware cleaning:
 - First with a mild detergent solution: a precaution for organic substances that are lipophilic
 - Final rinsing with DI water a couple of times:
 - 1st rinse with a very small volume
 - 2nd rinse with a somewhat bigger volume
- The best drying of water-rinsed containers and equipment is air drying. You can usually properly place inverted glassware and other equipment on the angled peg drying racks.