

*Chemistry 3A*

Introductory General Chemistry

*Experiment 7a*

Calorimetry

# Introduction

- **Calorimetry** means the measurement of heat. This measurement is of heat energy either given off (exothermic) or absorbed (endothermic) by the system from the surroundings
- Heat is measured by temperature changes with a thermometer and utilizes the specific heat capacity equation:


$$q = m \times c \times \Delta T$$

Symbol	Meaning	Units
<b><i>q</i></b>	energy as/in the form of heat	J or cal
<b><i>m</i></b>	mass of the substance in the system	g
<b><i>c</i></b>	specific heat capacity of substance: this is a physical property determined by scientists doing experiments	J / (g °C)
<b><i>ΔT</i></b>	change in temperature in °C: this is the final temperature minus the initial temperature	°C

# Background: Part 1

When a solid is dissolved in a liquid (solvent), there can be an energy change related to the molecules of the solid interacting with the molecules of the solvent.

- When heat is **released** in dissolving, reaction is **exothermic**
- When heat is **absorbed** in dissolving, reaction is **endothermic**






# Background: Part 2

Objects with different temperatures in contact with each other, heat from the warmer object transfers to the cooler object until both reach same temperature.

- Since energy is conserved, if  $\text{H}_2\text{O}$  as ice is put in contact with  $\text{H}_2\text{O}$  as water, it should be possible to measure temperature changes in the water as ice melts to see if this law is observed

# Equipment You Will Use



# Consumables



# Procedure: Part 1

1. Construct the calorimetry setup with the foam cups
2. Measure the mass of the cups ONLY
3. Add 45.0 to 55.0 mL DI water to grad cylinder & pour into cups
4. Record temperature to nearest 0.1°C
5. On glassine paper, scoop out 4.8-5.2 g citric acid without causing excessive waste. Tighten the reagent lid after.
6. Add citric acid carefully but quickly to the calorimeter with water and place lid on setup; VERY CAREFULLY with lid closed stir the solution until solid dissolves

# Procedure: Part 1

7. Record temperature to nearest 0.1°C after it stabilizes (it might go up or down)
8. Remove the cardboard lid. Record the mass of the calorimeter setup: cups + solution
9. Empty the solution into a beaker (for waste) then rinse and dry the inner cup
10. Switch the inner cup with outer cup
11. Repeat steps 3-8 with sodium carbonate ( $\text{Na}_2\text{CO}_3$ )
12. Pour out the waste into the beaker, discard the solution into the sink



# Procedure: Part 2

1. Construct the calorimetry setup with the foam cups and add 45.0-55.0 mL DI water
2. Measure the mass of calorimeter with water without lid
3. Record temperature to nearest 0.1°C
4. Transfer 2-3 cubes ice to calorimeter
5. Place lid and stir carefully with thermometer until ice melts. Record the final temperature
6. Record mass of calorimeter with contents (not lid)

# Example Data Analysis

## PART 1 – Exothermic and Endothermic Changes

### DATA

Mass of calorimeter (empty): 16.4831 g

	$\text{H}_3\text{C}_6\text{H}_5\text{O}_7$	$\text{Na}_2\text{CO}_3$
Initial Temperature ( $^{\circ}\text{C}$ )	22.8 $^{\circ}\text{C}$	22.1 $^{\circ}\text{C}$
Final Temperature ( $^{\circ}\text{C}$ )	19.3 $^{\circ}\text{C}$	27.9 $^{\circ}\text{C}$
Mass of calorimeter and solution (g)	72.0944 g	68.3891 g

## PART 2 – Specific Heat Capacity and Enthalpy of Fusion

### DATA

	Data
Mass of calorimeter + water (g)	64.3499 g
Initial Temperature of water ( $^{\circ}\text{C}$ )	21.4 $^{\circ}\text{C}$
Final Temperature of water ( $^{\circ}\text{C}$ )	15.8 $^{\circ}\text{C}$
Mass of calorimeter + water + ice (g)	67.8851 g

# Example Data Analysis

## PART 1 – Exothermic and Endothermic Changes

### DATA

Mass of calorimeter (empty): 16.4831 g

	H <sub>3</sub> C <sub>6</sub> H <sub>5</sub> O <sub>7</sub>	Na <sub>2</sub> CO <sub>3</sub>
Initial Temperature (°C)	22.8°C	22.1°C
Final Temperature (°C)	19.3°C	27.9°C
Mass of calorimeter and solution (g)	72.0944 g	68.3891 g

### CALCULATIONS – PART 1 – EXOTHERMIC AND ENDOTHERMIC CHANGES

\* Show your work, complete with units for H<sub>3</sub>C<sub>6</sub>H<sub>5</sub>O<sub>7</sub>. I will assume you did Na<sub>2</sub>CO<sub>3</sub> the same way

1. ☐ Mass of the solution

**How do you calculate?**

2. ☐ Change in Temperature of Solution

3. ☐ Heat of Solution (q<sub>solution</sub>)

### RESULTS – PART 1 – ENDOTHERMIC AND EXOTHERMIC CHANGES

	H <sub>3</sub> C <sub>6</sub> H <sub>5</sub> O <sub>7</sub>	Na <sub>2</sub> CO <sub>3</sub>
Mass of Solution (m)	<b>How do you calculate?</b>	
Specific Heat of Solution (c)	4.184 J/g°C	4.184 J/g°C
Change in Temperature of Solution(ΔT)		
Heat of solution (q <sub>solution</sub> )		
Heat of reaction* (q <sub>reaction</sub> + q <sub>solution</sub> = 0)		
Reaction is: Exothermic or Endothermic		

\* The water in the calorimeter is considered the reaction's surroundings, therefore if heat is released by the reaction, it is absorbed by the water and vice-versa.

# Example Data Analysis

## PART 2 – Specific Heat Capacity and Enthalpy of Fusion

### DATA

	Data
Mass of calorimeter + water (g)	64.3499 g
Initial Temperature of water (°C)	21.4°C
Final Temperature of water (°C)	15.8°C
Mass of calorimeter + water + ice (g)	67.8851 g

### CALCULATIONS – PART 2 – SPECIFIC HEAT CAPACITY AND ENTHALPY OF FUSION

\* Show your work, complete with units.

1. Mass of H<sub>2</sub>O

**How do you calculate?**

2. Change in Temperature of Water

**How do you calculate?**

3.  $q_{\text{water}}$

4. Mass of Ice

5. Moles of Ice

6.  $q_{\text{ice}}$

7.  $q_{\text{water}} + q_{\text{ice}}$

### RESULTS

	Results
Mass of H <sub>2</sub> O	
Specific Heat of water (c)	4.184 J/g°C
Change in Temperature of water ( $\Delta T$ ) ( $\Delta T = T_f - T_i$ )	
$q_{\text{water}}$ (kJ)	<b>How do you calculate?</b>
Mass of Ice	
Moles of Ice (n)	
Enthalpy of Fusion ( $\Delta H_{\text{fus}}$ )	6.01 kJ/mol
$q_{\text{ice}}$ (kJ)	
$q_{\text{water}} + q_{\text{ice}}$	

# Example Data Analysis

## PART 1 – Exothermic and Endothermic Changes

### DATA

Mass of calorimeter (empty): 16.4831 g

	H <sub>3</sub> C <sub>6</sub> H <sub>5</sub> O <sub>7</sub>	Na <sub>2</sub> CO <sub>3</sub>
Initial Temperature (°C)	22.8°C	22.1°C
Final Temperature (°C)	19.3°C	27.9°C
Mass of calorimeter and solution (g)	72.0944 g	68.3891 g

### CALCULATIONS – PART 1 – EXOTHERMIC AND ENDOTHERMIC CHANGES

\* Show your work, complete with units for H<sub>3</sub>C<sub>6</sub>H<sub>5</sub>O<sub>7</sub>. I will assume you did Na<sub>2</sub>CO<sub>3</sub> the same way

1. ☐ Mass of the solution 72.0944 - 16.4831 = 55.6553 g

4.184 J/g°C × 55.6553 g × (-3.5°C)  
= -820 J      2 significant digits

3. ☐ Heat of Solution (q<sub>solution</sub>)

2. ☐ Change in Temperature of Solution

19.3 - 22.8 = -3.5°C

### RESULTS – PART 1 – ENDOTHERMIC AND EXOTHERMIC CHANGES

	H <sub>3</sub> C <sub>6</sub> H <sub>5</sub> O <sub>7</sub>	Na <sub>2</sub> CO <sub>3</sub>
Mass of Solution (m)	55.6553 g	51.9060 g
Specific Heat of Solution (c)	4.184 J/g°C	4.184 J/g°C
Change in Temperature of Solution(ΔT)	-3.5°C	5.8°C
Heat of solution (q <sub>solution</sub> )	-820 J	1300 J
Heat of reaction* (q <sub>reaction</sub> + q <sub>solution</sub> = 0)	820 J	-1300 J
Reaction is: Exothermic or Endothermic	endothermic	exothermic

\* The water in the calorimeter is considered the reaction's surroundings, therefore if heat is released by the reaction, it is absorbed by the water and vice-versa.

# Example Data Analysis

## PART 2 – Specific Heat Capacity and Enthalpy of Fusion

### DATA

	Data
Mass of calorimeter + water (g)	64.3499 g
Initial Temperature of water (°C)	21.4°C
Final Temperature of water (°C)	15.8°C
Mass of calorimeter + water + ice (g)	67.8851 g

### CALCULATIONS – PART 2 – SPECIFIC HEAT CAPACITY AND ENTHALPY OF FUSION

\* Show your work, complete with units.

1. Mass of H<sub>2</sub>O

$$64.3499 - 16.4831 = 47.8668 \text{ g}$$

2. Change in Temperature of Water

$$15.8 - 21.4 = -5.6^{\circ}\text{C}$$

3.  $q_{\text{water}}$

$$\begin{aligned} 47.8668 \text{ g} \times (-5.6^{\circ}\text{C}) \times 4.184 \text{ J/g}^{\circ}\text{C} \\ = -1.1 \text{ kJ} \\ = 1.1215 \text{ kJ (no intermediate round off)} \end{aligned}$$

4. Mass of Ice

$$67.8851 - 64.3499 = 3.5352 \text{ g}$$

5. Moles of Ice

$$\begin{aligned} 3.5352 \text{ g} / 18.02 \text{ g/mol} \\ = 0.1962 \text{ mol} \end{aligned}$$

6.  $q_{\text{ice}}$

$$\begin{aligned} 0.1962 \text{ mol} \times 6.01 \text{ kJ/mol} \\ = 1.179 \text{ kJ} \end{aligned}$$

7.  $q_{\text{water}} + q_{\text{ice}}$

$$\begin{aligned} -1100 + 1179 &= 79 \text{ J} \\ -1121.5 + 1179 &= 57 \text{ J} \end{aligned}$$

# Example Data Analysis

## PART 2 – Specific Heat Capacity and Enthalpy of Fusion

### DATA

	Data
Mass of calorimeter + water (g)	64.3499 g
Initial Temperature of water (°C)	21.4°C
Final Temperature of water (°C)	15.8°C
Mass of calorimeter + water + ice (g)	67.8851 g

### RESULTS

	Results
Mass of H <sub>2</sub> O	47.8668 g
Specific Heat of water (c)	4.184 J/g°C
Change in Temperature of water ( $\Delta T$ ) ( $\Delta T = T_f - T_i$ )	-5.6°C
$q_{\text{water}}$ (kJ)	-1.122 kJ
Mass of Ice	3.5352 g
Moles of Ice (n)	0.1962 mol
Enthalpy of Fusion ( $\Delta H_{fus}$ )	6.01 kJ/mol
$q_{\text{ice}}$ (kJ)	1.179 kJ
$q_{\text{water}} + q_{\text{ice}}$	57 J

# Clean Up

- Return equipment cleaned if necessary to its storage areas
- Return cardboard lid to supply bin
- Discard the Styrofoam cups